

Key Features

- Create batches for each kind of antibody testing.
- Use any time period for the batch - day, week, month, or no time period.
- Assign batches to a technologist or a group of technologists so the system can still track testing assignment and results entry.
- Review test results as a batch to speed the review/sign out process.
- Select a no-review option and report results immediately after entry.
- Print a report of the batched tests or of individual results for a patient.

...and more

The **Antibody Screening for Large Volume Laboratories Module** makes batching all those tests for this month's antibody testing easy. Use batches to group tests by assay method, for a time period (if appropriate). Create a batch size that fits the method of testing you use.

Then review the results as a batch for speedy turn-around-times.

Use with Test Ordering by Business Rules to automatically order tests in batches.


1. Create a batch by applying a name, determining the batch size, and assigning it to an individual or group of technologists.
2. Move tests from the worklist into a batch.
3. Complete testing and select the batch for results entry.
4. Sign out immediately after entry or complete a one- or two-step review cycle.

System Requirements

- Microsoft Windows 95+
- Microsoft Access 97+ or SQLServer 7.0+
- 20 MB of disk space
- 32 MB or memory
- VGA (800X600)

For More Information

For more information about HistoTrac or for a free test-drive of HistoTrac, visit www.systemlink-inc.com. To contact a HistoTrac sales representative, call 540.548.2699, or send an e-mail message to histotrac@systemlink-inc.com



4230 Lafayette Center Dr.
Suite I
Chantilly, VA 20151
phone 703.295.0757
fax 540.785.3469
www.systemlink-inc.com