

### Key Features

- Print labels for tube mailers to obtain patient samples.
- Use barcodes to speed sample accessioning.
- Create mailing labels to the patient's home address or for batch mailers to a single address.
- Use printed labels for storing samples in inventory.
- Apply labels to tubes used during testing.  
...and more

The **Barcode and Label Modules** are available for 1) sending tube mailers to your patients, and 2) creating labels for test tubes used in the laboratory. Either module will save you time and reduce those errors. No more small print on small tubes with blunt-tipped markers! Every label is legible. And you chose the information you want on your labels.

If your phlebotomy department uses barcodes, speed high-volume sample accessioning by using a barcode reader.

Create barcode or regular tube labels and mailing labels for those sample tube mailers you send out for antibody screening samples.

### System Requirements

- Microsoft Windows 95+
- Microsoft Access 97+ or SQLServer 7.0+
- 20 MB of disk space
- 32 MB or memory
- VGA (800X600)

### For More Information

For more information about HistoTrac or for a free test-drive of HistoTrac, visit [www.systemlink-inc.com](http://www.systemlink-inc.com). To contact a HistoTrac sales representative, call 540.548.2699, or send an e-mail message to [histotrac@systemlink-inc.com](mailto:histotrac@systemlink-inc.com)

4230 Lafayette Center Dr.  
Suite I  
Chantilly, VA 20151  
phone 703.295.0757  
fax 540.785.3469  
[www.systemlink-inc.com](http://www.systemlink-inc.com)

