

Key Features

- Link members of the paternity case to view results together.
- Align haplotypes for easy viewing of inheritance patterns.
- Report the described haplotype assignments for easy interpretation by the legal representative.
- Store and report probability statistics.
- Keep the digital photo of each individual on file.
...and more

The **Paternity Testing Module** lets you document the required identification items needed to ensure legal presentation of the paternity case. Store HLA typing results. If you test red cell systems, use the Red Cell Typing Module alongside this module for complete test history.

1. Accession the samples on mom, child(ren), and alleged dad.
2. Collect driver's license #, Social Security #, race, date of birth, and a digital photograph to properly identify the individual.
3. Use the Worklist to track workflow.
4. Enter the HLA typing results of each individual and store them as a 'family'.
5. Print a Paternity Case Report.

System Requirements

- Microsoft Windows 95+
- Microsoft Access 97+ or SQLServer 7.0+
- 20 MB of disk space
- 32 MB or memory
- VGA (800X600)

For More Information

For more information about HistoTrac or for a free test-drive of HistoTrac, visit www.systemlink-inc.com. To contact a HistoTrac sales representative, call 540.548.2699, or send an e-mail message to histotrac@systemlink-inc.com

4230 Lafayette Center Dr.
Suite I
Chantilly, VA 20151
phone 703.295.0757
fax 540.785.3469
www.systemlink-inc.com

